

# MANUXELECT TRACTION



The traction version of the ManuXelect has a rolling section with a range of motion of 15 cm to perform traction therapy. This rolling section can also be locked in place for conventional use.


Your ManuXelect traction couch comes standard with a fixation rail on the head section and a traction frame for mounting the Eltrac 471 traction device.

The options for the height adjustment and the lifting capacity of the ManuXelect traction couch are the same as for the other ManuXelect versions.

Attention! The ManuXelect Traction is available only in the colours grey (colour number 013) and Sirius blue (colour number 004) and only with a section width of 67 cm.

## TRACTION ACCESSORIES

3xxx014	Flexistool grey, with separate height adjustment for each leg, available in 18 colours
3445022	Pelvic belt for Manumed Traction
3445021	Thoracic belt for Manumed Traction
1471808	Saunders clinical traction device including mounting bracket for the Eltrac 471
1671901	Eltrac 471, unit for traction therapy



## TECHNICAL SPECIFICATIONS

Size tabletop (lxw):	approximately 240 x 67 cm (incl. traction frame)
Movable leg section	15 cm
Lockable	yes

See brochure Eltrac 471 for more information.




# ACCESSORY: SAUNDERS CLINICAL TRACTION DEVICE

## CERVICAL TRACTION WITH THE SAUNDERS CLINICAL TRACTION DEVICE

Cervical traction can be very effective in the treatment of cervical hernias that require widening of the neural foramen or stretching of the soft tissue in the cervical structure. The Saunders clinical traction device makes cervical traction simple and comfortable.

With the Saunders clinical traction device, the traction forces are directed towards the occiput. The forces are not applied to the chin, so there is no pressure on the temporomandibular joint. The Saunders clinical traction device can be adjusted to individual needs so each patient receives the maximum benefit.

The angle of the cervical traction is adjustable with the Saunders clinical traction device. This permits optimal traction as described in the literature.



The combination of the Eltrac 471 with an Enraf-Nonius traction table and the Saunders clinical traction device is the ideal solution for anyone who is looking for responsible, safe and comfortable cervical traction.



### SAUNDERS CLINICAL TRACTION DEVICE

1471808 Saunders clinical traction device with mounting bracket

