

Mphi 5

the complete solution

The Mphi 5 has been designed to reach easily each area of the body and to convey the high quality energy of the MLS® pulse to a well-defined 5 cm diameter Target Area, whose shape and size has been optimized to guarantee an excellent therapeutic result. The size of the irradiated zone, in fact, makes it possible to treat homogeneously great tissue volumes, such as large fractions of muscle strips or whole joints. This causes a rapid response of the irradiated structures. The preset therapeutic protocols facilitate the operator to choose the kind of treatment, and the multidiode applicator Charlie offers the advantage to act over large areas of the body in short time, allows to work without requiring the constant assistance of an operator.

