

NEW POCKET PHYSIO IONOTENS



Description

The Pocket Physio Line includes a new electrostimulator for the most specific applications of pain relief and anti-inflammatory electrotherapy: the IONOTENS model. The device software contains 33 preset programs with the most used protocols plus 4 free programs for setting customized parameters; the dual power operation allows applications at home (with rechargeable battery) and in the therapy room (with power outlet). Due to its ease of use, personal use is possible for each patient.

It delivers antalgic tens currents, ionophoretic currents, anti-inflammatory and jontophoretic currents to relieve hyperhidrosis in hands and feet.

Basic equipment

- 1 set of pre-gelled electrodes (4 50x50 mm 2 50x90 mm)
- 2 fixing elastic bands
- 2 splitters
- 4 sponge bags (2 50x50 mm 2 60x80 mm)
- 4 conductive silicone electrodes (2 50x50 mm 2 60x80 mm)
- External power supply
- Rechargeable battery
- Output cable
- User manual and applications
- Suitcase

Programs / Applications list

- Cervicalgia
- Contractures
- Widespread bruises
- Sprains
- Sharp pain
- Chronic pain
- Gonalgia
- Iontophoresis
- Hyperhidrosis
- Metatarsalgia
- Sciatica
- Tore up
- Tendinitis

Technical features

Supply	rechargeable battery, external power supply
Wave forms	tens, iontophoresis
Number of output channels	2 independent
Current intensity (mA)	0-120
Frequency (Hz)	1-160
Issue	continuous / modulated
Pulse duration (us)	50-1s

PRODUCT CODE: NEWPF1201060

PRICE: € 650,00