



3D-adjustable and foldable side pelottes



Adjustable Therapy table side pelottes



3D-adjustable footplates



3D-adjustable knee supports



Novum size 3 with pelvis, thorax and knee supports and therapy table

## Indicated by the following foot malpositions:

- Weak or contracted foot malposition
- Equinus position and talipes calcaneus
- Talipes varus
- Other complicated foot malpositions

## Technical details: Novum

|  | Size 1       | Size 2       | Size 3        |
|--|--------------|--------------|---------------|
| Armpit height:                                     | 70 to 101 cm | 83 to 121 cm | 115 to 153 cm |
| Overall width:                                     | 59 cm        | 59 cm        | 70 cm         |
| Overall length:                                    | 85 cm        | 85 cm        | 112 cm        |
| Overall height*:                                   | max. 138 cm  | max. 153 cm  | max. 201 cm   |
| Overall height, horizontal:                        | 60 cm        | 70 cm        | 76 cm         |
| Standard support cushion – large (width x height): | 31 x 43 cm   | 36 x 50 cm   | 45 x 89 cm    |
| Standard support cushion – small (width x height): | 25 x 24 cm   | 32 x 26 cm   | 36 x 30 cm    |
| Tilt:  | -5 bis 90°   | -5 bis 90°   | -5 bis 90°    |
| Loading capacity:                                  | 50 kg        | 70 kg        | 100 kg        |

\*including all accessories (footplate to upper edge of supporting surface)