

SONO LE SOLUZIONI CREATIVE E TECNOLOGICHE
CHE FANNO UN PRODOTTO ORIGINALE



Maschera obiettivi



Maschera funzioni



Menù funzioni

Resistive and Capacitive Diathermy	STARTECAR 200 W
Power supply	115 - 230, 50/60 Hz, +/- 10%
Cabinet	Desk top cm 40x40x40 h
Coloured display	Yes with touch-screen 6,5"
Capacitive / Resistive	Yes, as standard
Continuous and pulsed emission modes	Yes
Working frequency	485 KHz +/- 10 Hz
Power max	200 W
Impedance on real-time	Yes
Temperature sensor	Yes
Hand-piece	Interactive Probe (or Ergonomic hand-piece on demand) *
Electrodes	Res/Cap, 30, 60 mm
Cream 1 Kg	Yes, as standard
Set of accessories for static applications	
60 mm electrode for Bipolar hand-piece	
Portable case	
Price with VAT	
Trolley mod. STAR	
Trolley mod. STAR Slim	

