

STAR TECAR



MAIN TECHNICAL FEATURES

- Power supply: 115– 230 Vac – 50/60 Hz
- Power consumption: max. 750 VA
- Protection Class: 1 BF
- EEC 93/42: IIa
- Therapeutic modalities: Resistive and Capacitive
- Application modalities: manual, bipolar, static
- Emission: thermal and non-thermal
- Working frequency: 480 KHz +/- 10%
- TOUCH control plus 10,1
- Impedance in real-time
- Temperature sensor for monitoring patient's feeling
- Pre-set programs
- Free memories
- Patients' data-base

Standard accessories: interactive Probe with LCD on board, n. 2 applicators (30 e 60 mm), bipolar hand-piece, 2 electrodes mod. FLEXI, neutral, elastic band, gell 250 gr

