

# Technical Specifications

## CPM FOR KNEE



1. Double Color touch screen control,dual display in English version
2. Knee sensor makes running angle more accurately.
3. Dual bars make it more stable ,which also have the function of reduction in total machine length .
4. Torque setting: 3 type
5. Double Screen setting of angle,speed,time
6. Special design of the support bar suiting to human curve
7. Input: 100 to 230 VAC , 50/60 Hz
8. Treatment Time: Continuous for 240 min
9. Speeds:0.5°-3.5°/s ( 9 levels)
10. Working Angles of knee joint: -5°- 120°
11. Suit to people 1.55m- 1.9m